


APPENDIX 1

400MHZ H1

Current Data Parameters
 NAME EXPNO 100
 PROCNO 1
 F2 - Acquisition Parameters
 Date_ 2002013
 Time 6.24
 INSTRUM spect
 PROBHD 5 mm QNP 1H
 PULPROG zgpg30
 TO 32768
 SOLVENT Aceton
 NS 32
 DS 0
 SWH 7163.508
 FIDRES 0.219235
 AQ 2.2607028
 RG 328.1
 DM 69.000
 DE 8.00
 TE 300.0
 D1 1.0000000
 D11
 ***** CHANNEL f1
 NUC1 1H
 P1 11.30
 PL1 0.00
 SFO1 400.132605
 F2 - Processing parameters
 SI 32768
 SF 400.1300535
 WDW EM
 SSB 0
 LB 0.00
 GB 0
 PC 1.00
 ID INEFT plot parameters
 CX 25.00
 CY 10.000
 F1P 4091.30
 F2P -0.500
 F3P -200.07
 SPACCH 0.42000
 MICH 168.05160

7.51738
 7.51334
 7.49664
 7.44177
 7.43808
 7.41998
 7.36645
 7.34774
 7.34039
 7.33589
 7.33213
 7.32306
 7.31887
 7.30386
 7.27365
 7.25566
 5.66520
 5.49553
 4.22873
 4.22244
 4.21105
 3.89657
 3.89283
 3.87487
 3.39203
 3.36688
 3.34896
 3.31376
 3.29602
 3.28814
 3.27040
 3.15445
 3.14309
 3.12881
 3.11748
 3.09821
 3.07626
 3.07325
 3.05122
 2.49540
 2.49088
 2.48640




APPENDIX 2

400MHz C13 (250p to -10p)

```

Current Data Parameters
NAME C1C-HBO-S1
EXPNO 400
PROCNO 2
F2 - Acquisition Parameters
Date_ 20030518
Time 13.27
INSTRUM spect
PROBHD 5 mm QNP 1H
PULPROG zgpg30
TD 65536
FIDRES 0.359443
AQ 1.2317875
RG 20842.5
DM 19.100
DE 10.00
TE 300.0
D1 2.00000000
d11 0.03000000
d12 0.00020000
***** CHANNEL f1
NUC1 13C
P1 9.50
PL1 -3.00
SFO1 100.6248023
***** CHANNEL f2
COPROG2 waltz16
NUC2 1H
PCPD2 88.00
PL2 120.00
PL12 10.00
PL13 22.100
SFO2 400.1318000
F2 - Processing parameters
SI 32768
SF 100.6127895
WDW EN
SSB 0
LB 1.00
GB 0
PC 0.60
10 MHz plot parameters
CX 25.00
FIP 180.000
F1 1810.30
F2 0.000
F2P 0.000
PPOH 7.20000
MICH 724.4187

```



S2/CDC13,1H(5mm)

Archive directory: /export/home/vnmr1/vnmrSYS/data
Sample directory:
File: PROTON

Pulse Sequence: s2pul

Solvent: CDCl3

Temp: 30.0 C / 303.1 K

INDVA-500 "BNEC500"

Relax. delay 1.500 sec

Pulse 45.0 degrees

Acq. time 3.000 sec

Width 8000.0 Hz

16 repetitions

OBSERVE H1, 499.8436750 MHz


DATA PROCESSING

Line broadening 0.1 Hz

FT size 65536

Total time 1 min, 12 sec

APPENDIX 3



S2/CDC13, 13C(5mm)

Pulse Sequence: s2pul

Solvent: CDC13

Temp. 30.0 C / 303.1 K

User: 1-14-87

INOVA-500 "BMECS00"

Relax. delay 2.000 sec

Pulse 42.4 degrees

Acq. time 1.042 sec

Width 3146.5 Hz

2000 repetitions

OBSERVE C13, 125.6857896 MHz

DECOUPLE H1, 499.8461695 MHz

Power 39 dB

Continuously on

WALTZ-16 modulated

DATA PROCESSING

Line broadening 1.0 Hz

FT size 131072

Total time 6 hr, 46 min, 47 sec

APPENDIX 4

33.428


42.995

65.268
66.009


76.748
77.256
77.980

126.417
126.150
126.417
126.574
126.685
126.792
126.939
138.939

158.681



APPENDIX 5



	RT	Area	% Area	Height
1	13.797	556442	100.00	23734

Analytical condition

HPLC Type:

Pump : Waters 600E
 Detector : Waters 2996 Photodiode Array Detector
 Autosampler: Waters 717 plus
 Mobile phase : 1% TEA, pH7.5 / MeOH = 40 / 60
 Flow rate : 1.0 mL/min
 Column : Inertsil 5 ODS-80A, 3.2*250-mm
 Column Oven : 40C
 Wavelength : 254 nm